

# STERIS 4085 GENERAL SUGICAL TABLE

ReVive Remanufactured by Bryton



**AMS-4085**  
General Surgical Table  
Remanufactured through ReVive  
by Bryton Corp.  
(Shown w/X-Ray Tops - Not included)

## AMS-4085 Specifications

GENERAL	CAPACITY
Operating Surface	81" L x 20"W (2057 x 508 cm)
Table Weight	820 lb (372 kg)
Max Permissible Load (Lift)	1100 lb (499 kg)
Max Permissible Load (Articulation)	600 lb (272 kg)

HAND CONTROLS	ADJUSTMENT
Height (w/o cushion)	26" - 45"W (660 x 1143 cm)
Lateral Tilt (left/right)	20°
Trendelenburg	30°±1°
Reverse Trendelenburg	30°±1°
Slide Motion	Head: 9"(227 mm) Foot: 9"(227mm) from center
Powered Kidney Bridge	4" (110mm)
Back Up/Down	40°/80°
Leg Up/Down	0°/105°
Flex/Reflex	140°/100°±1°
Base Braking (Lock/Unlock)	Yes
Return to Level	Yes

The STERIS 4085 General Surgical Table is a mobile, electro-hydraulically operated surgical table designed to support virtually all general surgical procedures including cardiac and vascular, endoscopic, gynecology, urology, nephrectomy, neurology, ophthalmologic and orthopedics with the addition of STERIS table accessories.

## Accessories & Pads

PRODUCT	OPTIONS & DESCRIPTION	QTY
Cushion Pad Sets for AMS-5085	<input type="checkbox"/> SC-3380 (Standard Pad) <input type="checkbox"/> PC-3380 (Premiere Pad) <input type="checkbox"/> IGC-3380 (Integra-Gel Pad)	<input type="text"/>
ACC-P2	Standard Accessory Package: 2 armboards w/2" pad (AS-4000), 2 armboard restraint straps (PR-4000), and 1 Standard patient restraint strap (PR-3200)	<input type="text"/>
ACC-CYSTO/URO	Cysto/Uro Accessory Package: 1 pair of Remanufactured Yellofins (LS-3000-SR), 2 Armboards w/ 2" Pads (AS-4000), 1 Patient Restraint Strap (PR-3200), 1 Drainage Bag System (DB-1000), and 1 Box of 10 Drainage Bags (DB-1055)	<input type="text"/>
ACC-URO	Uro Accessory Package: Set of Armboards, Drainage Bag System with 1 Box of 10 Drainage Bags.	<input type="text"/>

## Customer Info


First Name		Last Name	
E-Mail		Telephone	
Company		Cus. Account #	

### REVIVE WARRANTY

2 YEARS PARTS

2 YEARS LABOR

# AMS-4085 RANGE OF MOTION (TYPICAL)



Drawings are not to scale.  
Diagrams are to only illustrate  
typical range of motions.

<b>REVIVE WARRANTY</b>	
2 YEARS PARTS	2 YEARS LABOR